

S. Manimozhi

Founder & President Marketing

+91 9962514402 - 044 24954592

heartcare1995@gmail.com

www. heartcareindia.com



HEART LUNG MACHINE VERSION OFFER 1. INDIGENOUS



SPECIFICATION:

Heart Lung Machine Moveable Model 4 heads: feather touch with battery backup one hour indigenous and compact single designs ensure reliability conventional: electro mechanical drive continuous flow operation flow rate adjustable large & easy to read digital led display of flow in LPM & RPM hand crack for manual operation in case of power failure. Self – Locking control occasions three size foraging for use. In pump heat ½", 3/8", ¼" interchange ability of pump heats can be done quickly and easily , full stainless steel body timer, pressure gauge imported spot light(1yr warranty).



HEART LUNG MACHINE OFFER 2. INDIGENOUS WITH ALL STANDARD & OPTIONAL ACCESSORIES





	COMPLIANCE	
SPECIFICATIONS OF THE HEART LUNG MACHINE	YES/NO	
1. Heart lung is a apparatus, through which the blood is temporary		
diverted, during heart surgery, to oxygenate it and pump it	YES	
throughout the body to maintaining circulation until the heart lung	1 E3	
are able to return to normal function		
2. Operational requirements :	YES	
Basic equipment will consist of the following unit	YES	
a. 5 pump single head or 3 pump single head + 2 pump double		
head with base.	YES	
b. monitor	YES	
c. pressure monitor -arterial Cardioplega transducers	YES	
d . timer -at least 3 timer	YES	
e. Temperature monitor with at least 2 probes.	YES	
f. Cardioplega pump should have display of total volume of each		
infusion along with delivery time	YES	
g. Safety device	YES	



h. level sensor	YES
i. ultra sonic air center (optional)	YES
3. Accessories with include:	YES
a. Stainless steel line clamp 10 no's	YES
b . Stainless steel intra cardiac suckers 4 adult 4 pediatric	YES
4. Techinuical specification:	YES
1. Pump console	YES
2. The unit has five pump console compactly arranged with separate power	
supply and control modules. Have easy access connectors for interchanging the	
pump	YES
3. Each individual roller pump should be capable of running independently at	
available voltage	YES
4. It has spill proof base	YES
5. This unit will be supplied with a battery backup for at least two pumps, all	
safety systems and accessories will be functioning for 60 min. Switching over	
from power to battery backup will be automatically done. The battery unit will	
be built into the pump base and it will be recharged automatically when the system	
is operating with main power supply	YES



6. Our machine has individual pump head, has Harvey roller pumps with facility for tubing to be used adjustable from 1/4 " to 5/8" through 3/8" and 1/2 " by easily changeable mechanism	YES
7. Individual pump head has display in digital - the total infusion in liters and delivery time the flow rate in LPM and in RPM	YES
8. Each pump has easy mechanism for occlusions system for different thickness of tube available in the market , 1/32" to 3/32 "	YES
9. We provided hand crank facility as critical safety feature hand crank loading from top for fastest access.	YES
10. The console has a compact base mount for the entire pump head together, with pole and handle s	YES
11. We have variable, changeable tubing holder in each pump head :1/4", 3/8",1/2" and 5/8"	YES
12. We are providing removal oxygen holder	YES
13. Roller pump have a self diagnostic circuit with provision to dedicate and display critical alarm condition, Optional pulsatile mode with the	MEG
mounted on any of the blood pumps 14. It has monitor mount with adjustable monitoring arm.	YES YES



Monitor:	YES
1. Pressure monitor : facility to monitor one arterial line pressure	
and one Cardioplega monitor pressure (total 2)	YES
2. Along with necessary transducers, Cardioplega line pressure (total2); along with necessary with pressure transducers and also with accurate display and alarm facilities' audio and visual	
delivery with stop, reset, and start function	YES
Temperature :	YES
1. Four temperature displays for patient monitoring and for Cardioplega monitoring with digital display in Celsius should be	
available.	YES



Hypothermia Machine Offer 1. INDIGENOUS



SPECIFICATION:

HYPOTHERMIA MACHINE MOVEABLE MODEL: Feather touch indigenous and compact dual reservoir collar/heater fiber tank digital display. Hypothermia cooler/ Heater will develop and maintain temperatures between 2'C to 42'V, the unit is compact and occupies less flour space and moves on conductive swivel castors build circuit barkers enhances safety. Full stainless steel body.



HEATER COOLER MACHINE FOR CARDIOPULMONARY BY PASS NEW VERSION OFFER 2 BRAND HEART CARE. INDIGENOUS







HEATER COOLER MACHINE FOR CARDIOPULMONARY BY PASS NEW VERSION OFFER 2 BRAND HEART CARE. INDIGENOUS





SPECIFICATION:

- •Heater cooler device used during cardiothoracic surgeries to warm or cool a patient.
- •2 Separate closed water tanks-both having independent cooling and warming.
- •2 Separate heating coil with 2000 Watts capacity.
- •Independent temperature control of the patient and Cardioplega via separated water circuits
- •Silent operation noise level < 60 Db (A)
- •Gradient mode for physiologically optimized heating and cooling (OPTIONAL)
- •Polished stainless-steel cabinet easy to clean and damage, rust resistant.
- •Easy to maneuver using smooth running castor wheels with foot-lever operated brakes.
- •Temperature control range 2°C to 42°C
- •Flow capacity, patient water circuit 15-20 1/min.
- •Flow capacity, Cardioplega water circuit 9-11 1/min.
- •Minimum water pressure inside the heat exchanger 0.3-0.7 bar
- •Maximum water pressure inside the exchanger 0.3-0.7 bar.
- •Heating/cooling blanket (adult-1,pediatric-1)
- •Tubing set for heating/cooling blankets.
- •Adapter set for heating/cooling blankets.
- •No water shunt &spillage
- •Supply voltage 240 V ~ = 10% 60Hz, AC:50-60Hz



Temperature control New version Model (HYPOTHERMIA MACHINE):	COMPLIANCE YES/NO
1. Temperature control and monitor system with Cardioplega supply	YES
2. Simultaneous delivery of water arterial and Cardioplega heat exchanges and to thermal blankets to be available from suitable ports.	YES
3. To work with power supply 220 v 20 v 50 hz	YES
4.Pessure regulator blanket ports maintaining the temperature of arterial port	YES
5. temperature display rance 0 to 50 degree Celsius	YES
6. microprocessor based unit to control, cool, rewarm and maintain the temperature.	YES
7. water outlet temperature of heat exchanger and blanket range 0 -40 d c	YES
8. maximum flow performance oxygenator heat exchanger supply prot 15-22lpm for fast cooling, 480mmHg Maximum pressure, blanket 1.5 to 2.5 lpm at zero head.	YES
9. it is capable of providing ice water for Cardioplega independently with variable cooling rate \	YES
10. temperature probe model for the operating ranges 0-42 d c	YES
11. temperature probe to fit in standard oxygenator (bubble / membrane)	YES
12. this equipment has touch screen system with optional accessories of remote control model. Optional remits control unit should be capable of taking 3 temp. probes and display temp in digital readouts. Alarm limits setting for at least 3 probes at crucial sites.	YES



AIR OXYGEN BLENDER BRAND HEART CARE. INDIGENOUS Low-flow Blender With 11 l pm Cascade flow meter



SPECIFICATION:

Air oxygen blender for heart lung/ neonatal care has been specially designed with 10 lpm and 1lpm cascade flow meter at a resolution of 100ml. ± 3% FiO2 accuracy at any set point and oxygen concentration control 21 to 100%. It features gas supply failure alarm and has a minimum bleed of 2.5 lpm facilitated with bleed on/off facility.



AIR OXYGEN BLENDER BRAND HEART CARE. INDIGENOUS Features & accessories. Low-flow Blender With 0-11 l pm Cascade flow meter







Oxygen inlet hose & Air inlet Air pressure regulato Blender pole hole

hose



Oxygen pressure



AIR COMPRESSOR BRAND HEART CARE. INDIGENOUS



SPECIFICATION:

Air Compressor: Air compressor sucks the atmospheric air and delivery the oil free air minimum pressure of air is 50 pound sans maximum pressure of 60 pound most of important items in the air compressor like. Coolly coil fan safety air storage tank are made in medical stansred stainless steel. Safety filter air compressor outlet should connect to the gas blender.



MEDICAL AIR COMPRESSOR BRAND HEART CARE. INDIGENOUS



SPECIFICATION:

- Medical Air Source
- Compact
- Quick operation
- •Medical grade Oil free
- Micro filter dry air
- •High Flow Capabilities, Greater versatility



FIBRILLATOR WITH PEDAL BRAND HEART CARE. INDIGENOUS



High Realiability and Safety
Adjustable Fibrillation Current
Light and Compact
Input 230 V and Output 0-12V
The electronic Fibrillator is specially
designed for Open Heart
Surgical Procedures
The cardiac Fibrillator provides and
adjustable alternating current to fibrillate
the heart.

The patient electode is of the 2-lead isolated type and can be easily autoclaved/sterlised. The fibrillation current delivered to the patient is controlled by a rotary potentiometer and displayed on a panel meter. The current



FIBRILLATOR WITH PEDAL BRAND HEART CARE. INDIGENOUS



Fibrillator With Pedal Specification ON/OFF SWITCH:

The switch serves to turn the unit on or off. The indicator in the switch lights up when the unit is turned on INPUT VOLTAGE 220V-1AMPS.

OUTPUT VOLTAGE POTENTIO METER

The output voltage is adjustable from 0-12V using Potentio meter should be set to OV when turning on the unit the voltage can then immediately be adjusted up to 12V. However, if the potentio meter had not been set to OV prior to turning on the unit the rest LED will start blinking. The unit is then ready for adjusting the output voltage.

OUTPUT CONNECTIONS FOR FIBRILLATOR ELECTRODE CONNECTORS:

The cable of the electrode has a length of approximately 3M banana Plug which are to be plugged into output sockets.



FIBRILLATOR WITH PEDAL BRAND HEART CARE. INDIGENOUS



POWER CABLE:

The power cable is firmly attached to the fibrillator has a length of 3M

PREPARING FOR OPERATION:

There are only few preparatory steps to take place the unit outside the sterile Area. The cable and the electrodes have been sterilized (if applicable for the Operation at hand) and ready for use.

- Check if the mains voltage complies with the voltage specified on the rear panel.
- Connect the power cable to a wall outet.
- Prior to operation, check safety circuit with power switch turned off, turn the Potentiometer to center position, turn power switch on the indicator.

OPERATING THE FIBRILLATOR:

- Ensure that potential meter is pointing to OV
- Turn power switch on the indicator.
- Slowly turn potential meter to increase output voltage clockwise until cardiac arrest is observed

Caution: Do not paddle joining together after use...keep it separately.



OPTIONAL ACCESSORY STAINLESS STEEL PEDAL OR ACRYLIC PEDAL BRAND HEART CARE. INDIGENOUS













REUSABLE CARDIOPLEGA DELIVERY SYSTEM ADULT & PEDIATRIC BRAND HEART CARE. INDIGENOUS



SPECIFICATION:

Helps to store to cooling the blood at the time of surgery, it can be reusable many no of times to multiple patient before the surgery the equipments should be taken care precaution and sterilized.


Note:

Comparing to other brand of Disposable Cardioplega the purchase price will be more than our brand



Reusable Cardio Plegia Cool Chamber and Silicon Washer Set. Optional Accessory







OUR BRANDS VS OTHER IMPORTED BRANDS:

	HEART LUNG MACHINE/HYPOTHERMIA MACHINE				
Heart Care Brands			Imported Brands		
•	We have a battery backup up to 2 hours for all 4 heads.	•	45 min of battery backup only avaliable		
•	Service and Maintaience, Spare parts are avaliable immediately with cost effective	•	Service and Maintaience, Spare Parts are not avaliable immediately with high cost along with duty charges		
•	We offer 3 years of spares warranty and 7 years Amc. Totally no of warranty is 10 years.	•	Fixed Total no of warranty maximum will be from 3 years without spare parts warranty		
•	Our Heads of the machine can be easily interchangeable	•	Heads are unable to interchange		
•	Once you fix the occulation you need not keep changing the settings on monthly basis for different patients without any interruption	•	Once you fix the occulation you need to keep setting up different features for on monthly basis for patients		
•	All price is very competative with all above feature	•	All price is not competative compare with our brand		
•	New version of Hypotherima you can do simutaneously warming the patient and cooling the heart, with touch screen and remote control	•	You Do not have warming the patient and cooling the heart simutaneously in any imported brand machine		
•	It works with our enviromental humidity temperature and on voltage fluctuation	•	It is very sensitive and non enviromental for our humidity temperature and voltagefluctuation especially in rual areas		
•	Any additional accessories like Blood sensor level, bubble dector, pressure monitor with probe are given in a affordable price.	•	Any additional accessories are available at a high cost compare without brand		



BLANKET FOR HYPOTHERMIA MACHINE WITH SUITABLE CONNECTORS (BRAND-CINCINNATI SUBZERO-U.S.A) (FDA - CE APPROVED)





DISPOSABLE BLANKETS/PADS

ADULT SIZE- 24*60 (61 CM*152.4CM)
PEDIATRIC SIZE- 24*60 (61 CM*152.4CM)
INFANT SIZE- 24*60 (61 CM*152.4CM)

RESUABLE/HYPOTHERIMA BLANKETS

ADULT SIZE- 24*60 (61 CM*152.4CM)
PEDIATRIC SIZE- 24*60 (61 CM*152.4CM)
INFANT SIZE- 24*60 (61 CM*152.4CM)



STERNUM SAW SET: Brand HCIC/BOJIN CE&ISO APPROVED





Specification:

Auto clavable
Reciprocating blade type
Ideal for sternotomy
Available in pediatric and
adult type

Contains : Hand piece, Saw guard, battery with charge , Aspetic channel, Sterlization Container Box, Saw Blade (5Nos). Warranty 1yr and Battery 6 Mnths.



REDO SAW SET: Brand HCIC/BOJIN



Specification:

Auto clavable
Reciprocating blade type
Ideal for sternotomy
Available in pediatric
and adult type

Contains : Hand piece, Saw guard, battery with charge , Aseptic channel, Sterilization Container Box, Oscillating Saw Blade (5Nos). Warranty 1yr and Battery 6 Months.



Head Band Univet- Brand Italy (CE approved)



Specification

Flexible structure with no pressure points adjustment of the vertical axis and declination angle interchangeable sweatband sponges can be worn over prescription eye wear, to suit with ESO HP headlight.



Head Light Univet- Brand Italy EOS HP LED Head Light



Specification:

Colour temperature 5700k
Brightness 45000lux
Battery life 8hrs max high density
Weight 28 g



Power Back System Univet- Brand Italy



Power back system for EOS HP LED HEADLIGHT without remote control CONTROL



Clip Univet- Brand Italy



Note:

Magnetic Clip to attach univet headlight any brand of loupe expect karl zeiss loupe



Headlight brand HCIC/CHAMP FOND CE APPROVED HEADLIGHT ONLY



Specification:

- Lamp weight is 12.9g
- •Light intensity 30,000lux@34cm, 21000lux@42cm
- •Spot size dia 9cm@42cm
- •Color temperature 5700k
- •Brightness level 10 levels
- Battery life 6 hrs highest intensity
- •Charging time 3-4 hours
- •Battery type rechargeable lithium ion
- •Battery capacity 600mah
- Power back weight 180g



HEADLIGHT WITH SUTIABLE HEAD BAND HCIC/CHAMP FOND CE APPROVED (the headlight attached the headband)



Note:

• Adjustable headband in all level



Clip BRAND HCIC/CHAMPFOND- CE APPROVED



Note:

 Adaptors clip can be attached any type of loupe except karl zeiss loupe

THANK YOU