



# MASTER BOOK

## 2023



### **HEART CARE INDIA CO**

No 21/11, 2<sup>nd</sup> Floor, VC Garden 2<sup>nd</sup> Street,  
Mandaveli, Chennai-600 028.

Phone No: 044-24954590/9962514402

Mail Id: [heartcare1995@gmail.com](mailto:heartcare1995@gmail.com)

Web: [www.heartcareindia.com](http://www.heartcareindia.com)

**Cardiac  
Vascular  
Cardiology**

# CARDIAC SURGERY

## CHASE HEART STABILIZER / SVR




## CHASE - INTRACORONARY VESSEL SHUNT


Bulb Tip short Shunt 30mm




Bulb Tip LONG Shunt 12mm



Tapered Tip Shunt 12mm



## TRIPOD



## Vessel Loops



## Pruitt F3® Balloon




## Flexcel™ 8-14 Fr



## XenoSure® Biologic Patch




## NIKO Flexible Pen




# CARDIAC SURGERY



LifeSpan ePTFE Graft



AlboGraft Polyester




OMNIFLOW



Aortic Punches

OTW – Embolectomy



Soft tissue Retractor

**BULLDOG CLAMP**



Wound Protector





# Viper II Vacuum Stabilization System



*Complete Line to Beating Heart Stabilization Technology*

*The Viper II Vacuum Stabilization System combines vacuum stabilization with an apical positioner and blower mister visualization.*



**BLM-700MS Viper II Blower Mister**

## Viper II Features

- *Low Profile Stabilization for Access to Anastomotic Site*
- *Malleable Footpads provide Optimal Positioning*
- *Apical Positioning for Access to All Arteries*
- *Blower Mister Offers Clear View and Access*
- *Stabilizer and Apical Positioner Fit Other Sternal Retractors*



## **HEART CARE INDIA CO**

**No 21/11, 2<sup>nd</sup> Floor, VC Garden 2<sup>nd</sup> Street,  
Mandaveli, Chennai-600 028.**

**Phone No: 044-24954590/9962514402**

**Mail Id: [heartcare1995@gmail.com](mailto:heartcare1995@gmail.com)**

**Web: [www.heartcareindia.com](http://www.heartcareindia.com)**



**MEDTRONIC equivalent to OCTOPUS 3 Stabilizer and URCHIN Heart Positioner**

**STB-700VS: Viper II Vacuum Stabilizer, M base**  
**API-700VS: Viper II Apical Lift Positioner, M Base**

**Designed to fit on the Medtronic or Chase Retractor**



**Vacuum Foot Pods and Vacuum Apical Positioner are identical on both "M" and "Q" bases. They are more equivalent to the Maquet system.**



**Disposable Retractor Rails to fit on the Maquet Retractor Sold as part of the STB-800VS: Viper II Vacuum Stabilizer with Q Base**

**New Codes: MAQUET EQUIVALENT**

**STB-800VS, Viper II Stabilizer with Q Base, includes Retractor Rails (see next image)**

**API-800VS, Viper II Apical Lift Positioner with Q Base**

**Equivalent to ACROBAT SUV, MODEL: OM-9000 Stabilizer and EXPOSE 4 Apical Positioner**

**Fits only on the Maquet Retractor**




# CHASE INTRACORONARY LONG 30 MM SHUNTS




**HEART CARE INDIA CO**

No 21/11, 2<sup>nd</sup> Floor, VC Garden 2<sup>nd</sup> Street,  
Mandaveli, Chennai-600 028.  
Phone No: 044-24954590/9962514402  
Mail Id: [heartcare1995@gmail.com](mailto:heartcare1995@gmail.com)  
Web: [www.heartcareindia.com](http://www.heartcareindia.com)



Bulb Tip  
LONG Shunt  
30mm



Bulb Tip Short  
Shunt  
12mm



Tapered Tip  
Short Shunt  
12mm

## TAPERED Tips

- Helps prevent Arrythmias
- Control arterial bleeding
- Faster & more accurate
- Soft Silicone Tip

## BULB Tips

- Minimize Trauma
- Obstruction –Free flow Path
- Improve Blood Flow
- Flexible Shaft

Product	Ref Code
Tapered Tip Shunts 1.00/ 12mm Shaft	SHT-1210S
Tapered Tip Shunts 1.00/ 12mm Shaft	SHT-1213S
Tapered Tip Shunts 1.00/ 12mm Shaft	SHT1215S
Tapered Tip Shunts 1.00/ 12mm Shaft	SHT-1218S
Tapered Tip Shunts 1.00/ 12mm Shaft	SHT-1220S
Tapered Tip Shunts 1.00/ 12mm Shaft	SHT-1225S
Bulb Tip Shunts 1.50/ 12mm Shaft	SHB-1215S
Bulb Tip Shunts 1.75/ 12mm Shaft	SHB-1218S
Bulb Tip Shunts 2.00/ 12mm Shaft	SHB-1220S
Bulb Tip Shunts 2.50/ 12mm Shaft	SHB-1225S
Bulb Tip Shunts 3.00/ 12mm Shaft	SHB-1230S
Bulb Tip Shunts 2.00 / 30mm Shaft	SHB-3020S
Bulb Tip Shunts 2.50 / 30mm Shaft	SHB-3025S
Bulb Tip Shunts 3.00 / 30mm Shaft	SHB-3030S

# LeMaitre- Range of Vascular Grafts



## LifeSpan® ePTFE Vascular Grafts

- Superb suture retention<sup>2</sup>
- High burst strength<sup>2</sup>
- for Haemodialysis



## AlboGraft® Polyester Vascular Graft

- Knitted and woven double velour grafts
- Approved collagen impregnation process
- Comfortable sealing
- Excellent suture ability
- Softness and conformability
- Removable external spiral for peripheral grafts



## Omniflow® II Biosynthetic Vascular Prosthesis

- Made Up of collagen and a polyester mesh  
For Hemodialysis where vein is not available



# XenoSure<sup>®</sup>

## Biologic Patch



**Exceptionally Strong, Uniform  
and Easy to Suture**

- Vascular & Cardiac Repair
- Suture Line Reinforcement
- Soft Tissue Deficiency Repair
- **NEW** – Dura Repair

**STRONG, UNIFORM AND  
EASY TO SUTURE**

- Biocompatible
- Excellent tensile and suture retention strength
- Does not require special sutures
- Uniform collagen thickness\* for easier suturing



*450K+ XenoSure Patches  
Implanted over 5 Years\*\**

\* Data on file at LeMaitre Vascular, Inc.

\*\* Based on LeMaitre Patch units sold in the last 5 years



**HEART CARE INDIA CO**

No 21/11, 2<sup>nd</sup> Floor, VC Garden 2<sup>nd</sup> Street,  
Mandaveli, Chennai-600 028.


Phone No: 044-24954590/9962514402

Mail Id: [heartcare1995@gmail.com](mailto:heartcare1995@gmail.com)


Web: [www.heartcareindia.com](http://www.heartcareindia.com)

# The SVR System *Left Ventricular Repair System*

The SVR<sub>v</sub> procedure is a comprehensive surgical therapy for patients suffering from ischemic heart failure. Studies have shown that surgically restoring a dilated left ventricle to a more normal elliptical anatomy using the SVR system results in significantly improved ejection fraction and less risk of post-operative mitral regurgitation



**Spherical Ventricle Post SVR**  
Undesirable surgical outcome



An SVR technique should result in a normal elliptical ventricle every time.

## Key Benefits:


- Standardizes the SVR procedure
- Achieves optimal ventricular function and volume
- Defines the new apex
- Reorients the mitral apparatus and papillary muscles properly



# Standardizing Surgical Ventricular Restoration (SVR) for consistent and elliptical results every time

**The SVR System includes:**

- The Mannequin™
- Defines the new apex
- Acts as a template to create an elliptical ventricle with optimal volume every time
- Offered in various sizes based on the patient's BSA
- Patch Seizer
  - Reduces potential for incorrectly trimming the patch
  - Decreases the procedure time




## Ordering Information

Product	Code	Description	Units
	SVC-0080S	SVR System with 80cc Mannequin, without patch	1
	SVC-0085S	SVR System with 85cc Mannequin, without patch	1
	SVC-0090S	SVR System with 90cc Mannequin, without patch	1
	SVC-0100S	SVR System with 100cc Mannequin, without patch	1
	SVC-0110S	SVR System with 110cc Mannequin, without patch	1
	SVC-0120S	SVR System with 120cc Mannequin, without patch	1
	SVC-0130S	SVR System with 130cc Mannequin, without patch	1
	SVC-0140S	SVR System with 140cc Mannequin, without patch	1
	SVC-0150S	SVR System with 150cc Mannequin, without patch	1

\* Patents US6681773 and additional patents pending

# CHASE -TRIPOD™

## Beating Heart Stabilization System



### *Easy to Use and Clinically Flexible:*

The TRIPOD™ Beating Heart Stabilization System combines ease of use with a rock solid, metal articulating arm and ultra-low profile foot pads. The strength and flexibility of the metal articulating arm delivers superior access and stabilization even to the most challenging Anastomotic locations.


### **TRIPOD Features:**

- Ease of use
- Superior arm strength and stability
- Low profile design
- Optional Atraumatic vessel occlusion



# Pruitt F3<sup>®</sup> Carotid Shunts

Cerebral Perfusion with Dual Balloon Occlusion



## Pruitt F3-S Polyurethane Outlying Carotid Shunt

Description	Diameter	Length	Model #
Shunt with T-Port	8F	31 cm	2015-10
Shunt with T-Port	9F	31 cm	2014-10

**Soft, flexible polyurethane with atraumatic tips to minimize intimal trauma.**

**Increased Flow: Average flow is 39% - 104%\* higher than leading competitors.**

**Improved Markings: Centimetres markings span the length of the shunt with a bar denoting centre point for easier identification of placement in carotid artery. Tether: Aids in removal of shunt.**

***\*Flow comparison testing conducted at body temperature at constant pressure..***

# Flexcel™ Carotid Shunt

Increased Flexibility and Flow

- Atraumatic tips
- Pre-attached tether
- Flexible material



## Flexcel Carotid Shunts (Inlying)


### Single Unit Packing

	Length	Model # <span style="border: 1px solid black; padding: 2px;">REF</span>
8F - packing with 5x1 unit 8F shunt	14.5 cm	2020-15
10F - packing with 5x1 unit 10F shunt	14.5 cm	2020-25
12F - packing with 5x1 unit 12F shunt	14.5 cm	2020-35
14F - packing with 5x1 unit 14F shunt	14.5 cm	2020-45



# NIKO AORTIC PUNCHES

## Single Rotating Head









### Aortic Punches:


1. Punches with Precision
2. Comfortable Grip
3. Smooth and easy to operate
4. Economical
5. Higher number of Usages

### Sizes:

- 2.8mm
- 3.6 mm
- 4.0 mm
- 4.4 mm

# NIKO POLYMER CLIP & MICRO CLIPS WITH OPEN CLIP APPLICATOR

	PLTX220	Medium	Blue ●
	PLTX230	Medium-Large	Green ●
	PLTX240	Large	Purple ●
	PLTX250	Extra-Large	Gold ●
	TCTL1201-24	CardioSet-ST	Red
	TCTL2200-24	CardioSet-MT	Blue



# SOFT TISSUE RETRACTOR

## Minimally Invasive Surgery

CARDIAC SURGERY / CTVS SURGERY / LAPROSCOPIC-GENERAL SURGERY  
GYNAECOLOGY / C-SECTION AND BREAST SURGERY



1. Assisted thoracic small incision surgery
2. Pediatric Congenital heart repair surgery
3. Minimally Invasive Cardio Surgery, valve replacemer
4. The small incision robot arm collaborative surgery



## HEART CARE INDIA CO

No 21/11, 2<sup>nd</sup> Floor, VC Garden 2<sup>nd</sup> Street,  
Mandaveli, Chennai-600 028.  
Phone No: 044-24954590/9962514402  
Mail Id: [heartcare1995@gmail.com](mailto:heartcare1995@gmail.com)  
Web: [www.heartcareindia.com](http://www.heartcareindia.com)



# BULLDOG CLAMP

**Vascu®Stop**



**Standard Silicone Coated  
Soft Coated Gentle Jaws**

**BSFS-009**



**BSFC-008**




**Multi Clamping Pressures**

## VASCUSTOP BULLDOG CLAMPS

Product Code	Products Description/ürün	barcode no	Pressure Gr.	In Box
BMXC-006	Maxi Curved	8680041680106	125/175/200 gr	25
BMXS-007	Maxi Streight	8680041680113	125/175/200 gr	25
BSFC- 008	Standard Curved	8680041682766	125/175/200 gr	25
BSFS-009	Standard Streight	8680041682773	125/175/200 gr	25
BMNC-010	Mini Curved	8680041682452	125/150 gr	25
BMNS-011	Mini Streight	8680041682759	125/150 gr	25
BSTC-029	Std.Soft Coated Curved	8680041680632	175/250 gr	25
BSTS-030	Std.Soft Coated Straight	8680041680625	175/250 gr	25

**Maxi,Standard & Mini Types**



# SILICON VESSEL TIES / VESSEL LOOPS

## BATRIK MEDICAL



**BLUE**

**RED**

**WHITE**

**YELLOW**



**Size**

**Mini ( 0.2cm \* 40cm)**

**Maxi (0.3 cm \* 40cm)**

**Medishare Inc,**

S-439, G.K.1, New Delhi ,

+91-11-29231218

Email-medishare@gmail.com

**Vessel loops are used in occluding, retracting and identifying arteries, veins, tendons and nerves in surgical procedures. Our vessels loops are extruded from medical grade silicone, with the exception of the Richard – Allan vessel loops, which are made of high performance, medical grade elastomeric. All vessels loops are latex free and radiopaque. Size Mini ( 0.2cm \* 40 cm & Maxi (0.3 cm \* 40cm)**



# The New Innovative ((70mm-200mm) NIKO Flexible Cautery with LED/USB

## Indications :

1. Onco Surgery / Abdominal Whipple
2. Liver Cancer Colorectal
3. Upper Spine
4. Hip Replacement / Hip Resection
5. Pelvic Resection / Hysterectomy
6. MICS, CABG, Neuro Endoscopy, LIMA
7. Kidney Transplant Surgery

**Flexible  
Pen**



## Features:

1. Easy to Operate
2. USB Connection with Long Cable
3. Ergonomically Design
4. Light in Weight
5. Adjust to any length
6. Long Blade ( 70mm-200mm)
7. Power Source 220ac/ Power Banks/ Laptop or any USB outlet
8. On the Spot Illumination

**Flexible  
Pen  
with light**






# Defibrillation Electrode

**Medishare Inc,**  
 S-439, G.K.1, New Delhi ,  
 medishare@gmail.com

303A5030		Connector compatible with: Phillips Heartstart: FR2, FR2 (+), FR3, Heartstream Forerunner, Heartstart MRx, Heartstart XL, Heartstart XL+
303A5035		Connector compatible with: Medtronic Physiocontrol, Osatu Bexen, Cardioline, Mindray, JOUSING
303A4530		Connector compatible with: Nihon Kohden:
303A2545		Connector compatible with: Zoll Medical Corp
303A6030		Connector compatible with: Zoll: E Series, M Series, R Series, AED PLUS, AED PRO
303A1525		Connector compatible with: CU Medical Systems, Inc, Cmos Drake, Progetti
303A2060		Connector compatible with: <b>Metrax Primedic</b>
303A4531		Connector compatible with: Nihon Kohden
303A7035		Connector compatible with: CU Medical Systems, Inc, Philips
303A4525		Connector compatible with: Schiller
303A4035		Connector compatible with: Mediana





# Niko

Electrosurgical Cautery Pad / Diathermy Pad /With & without Cables Valleylab / Bertchold / Erbe Connectors





STANDARD & REM



BUTTERFLY SHAPE



INFANT & NEONATAL



WITHOUT CABLE PAD

## CONNECTORS



1. Plates are available in two types SPLIT & STANDARD , These are hydro gel combined, very cool & comfortable
2. These are maximum temperature 4°C have no burning or marks on the skin
3. These are adaptable resistance.
4. These are available for Infant, Pediatric & Adult, with & without cable
5. These are most economical & use to avoid shock to the patient during surgery.




# NIKO Range of Electrosurgical Pencils

**NIKO PEN**




**NIKO PEN WITH NON STICK BLADE**



**NIKO PRO PENCIL WITH CERAMIC BLADE**



**NIKO REUSABLE PENCIL**



**NIKO FLEXIBLE PEN**



**NIKO FLEXIBLE WITH LIGHT**



**NIKO SMOKE PENCIL**



**NIKO PENCIL WITH TIP CLEANER**






# Vascular Surgery

# Pruitt F3<sup>®</sup> Carotid Shunts

Cerebral Perfusion with Dual Balloon Occlusion



## Pruitt F3-S Polyurethane Outlying Carotid Shunt

Description	Diameter	Length	Model #
Shunt with T-Port	8F	31 cm	2015-10
Shunt with T-Port	9F	31 cm	2014-10

**Soft, flexible polyurethane with atraumatic tips to minimize intimal trauma.**

**Increased Flow: Average flow is 39% - 104%\* higher than leading competitors.**

**Improved Markings: Centimetres markings span the length of the shunt with a bar denoting centre point for easier identification of placement in carotid artery. Tether: Aids in removal of shunt.**

***\*Flow comparison testing conducted at body temperature at constant pressure..***

# Flexcel™ Carotid Shunt

Increased Flexibility and Flow

- Atraumatic tips
- Pre-attached tether
- Flexible material




## Flexcel Carotid Shunts (Inlying)

### Single Unit Packing

	Length	Model # <span style="border: 1px solid black; padding: 0 2px;">REF</span>
8F - packing with 5x1 unit 8F shunt	14.5 cm	2020-15
10F - packing with 5x1 unit 10F shunt	14.5 cm	2020-25
12F - packing with 5x1 unit 12F shunt	14.5 cm	2020-35
14F - packing with 5x1 unit 14F shunt	14.5 cm	2020-45



# LeMaitre® Clot Management Catheters



Over-the-Wire Latex Embolectomy Catheter

Single Lumen Latex Embolectomy Catheter

NovaSil Single Lumen Silicone Embolectomy Catheter

**Over-the-wire Latex Embolectomy Catheter**

**Single Lumen Latex Embolectomy Catheter**

**Novasil Single Lumen Catheter**

**Latis Graft Cleaning Catheter**

Unique Latis mesh allows for efficient cleaning and effective removal of clots from synthetic grafts. The traction balloon engages tough adherent thrombi and maintains tactile feedback for effective treatment



**Syntel Catheter**

Shaft design ensures stability in center of vessel to minimize potential damage to the vessel intima. Non-fragmenting balloon enhances clinical outcomes.



## LeMaitre Embolectomy Catheters -

- Concentric balloons with progressive Inflation
- Radiopaque catheter shaft
- French size colour coding • Depth markings each 10 cm
- Robust Navigation

# XenoSure<sup>®</sup>

## Biologic Patch



Exceptionally Strong, Uniform  
and Easy to Suture

- Vascular & Cardiac Repair
- Suture Line Reinforcement
- Soft Tissue Deficiency Repair
- NEW – Dura Repair

**STRONG, UNIFORM AND  
EASY TO SUTURE**

- Biocompatible
- Excellent tensile and suture retention strength
- Does not require special sutures
- Uniform collagen thickness\* for easier suturing



450K+ XenoSure Patches  
Implanted over 5 Years\*\*



\* Data on file at LeMaitre Vascular, Inc.

\*\* Based on LeMaitre Patch units sold in the last 5 years



**HEART CARE INDIA CO**

No 21/11, 2<sup>nd</sup> Floor, VC Garden 2<sup>nd</sup> Street,  
Mandaveli, Chennai-600 028.  
Phone No: 044-24954590/9962514402  
Mail Id: [heartcare1995@gmail.com](mailto:heartcare1995@gmail.com)  
Web: [www.heartcareindia.com](http://www.heartcareindia.com)

# LeMaitre- Range of Vascular Grafts



## LifeSpan® ePTFE Vascular Grafts

- Superb suture retention<sup>2</sup>
- High burst strength<sup>2</sup>
- for Haemodialysis



## AlboGraft® Polyester Vascular Graft

- Knitted and woven double velour grafts
- Approved collagen impregnation process
- Comfortable sealing
- Excellent suture ability
- Softness and conformability
- Removable external spiral for peripheral grafts



## Omniflow® II Biosynthetic Vascular Prosthesis

- Made Up of collagen and a polyester mesh
- For Hemodialysis where vein is not available



# VenaBLOCK -Varicose Vein Closure System

VenaBLOCK® offers a unique treatment for venous insufficiency. Treatment is simply delivering cyanoacrylate inside of the dysfunctional vein segment. (VSM, VSP, or perforator veins) continuously while applying pressure over target vein segments



Catheter Length: **100 cm (PTFE)**

Catheter Material : Marked on every 2 cm  
OD: 2mm

Catheter Diameter **15 cm**

Catheter Shaft : **Braided, laser Guided**

Marker : Double Lumen , ID 0.5mm for agent ,  
ID 0.9mm for laser


# TemREN - Atherectomy System

Ready to use set content provides a practical application option. Rotation speed can be adjusted easily and provide safe space to surgeon.

## Consistent Solution For Peripheral Arterial Disease

Ready to use set content provides a practical application option. No additional installation equipment required before the process.



Rotation speed can be adjusted easily and provide safe space to surgeon.



## HEART CARE INDIA CO

No 21/11, 2<sup>nd</sup> Floor, VC Garden 2<sup>nd</sup> Street,  
Mandaveli, Chennai-600 028.  
Phone No: 044-24954590/9962514402  
Mail Id: [heartcare1995@gmail.com](mailto:heartcare1995@gmail.com)  
Web: [www.heartcareindia.com](http://www.heartcareindia.com)

# AnastoClip GC For Interrupted Vascular Anastomosis - AV Fistulas



## Key Benefits

- Interrupted, compliant anastomoticline
- Securely grips into vessel walls
- Fewer revisions to maintain functional patency<sup>1</sup>
- Reduced procedure time<sup>2</sup>

## Applications :

- AV fistulas
- Peripheral bypass vein grafting
- Arteriotomy and venotomy
- Vessel repair

Size	Code	Clip Width	Clip Height	Units per pack	Clip Length	Clip Width	Units per pack	Clip Length	Code
Medium	M	2.3 mm	1.1 mm	1 unit per pack	1.1 mm	2.3 mm	35	8 cm	4009-06
Medium	M	1.1 mm	1.1 mm	1 unit per pack	1.1 mm	2.3 mm	35	NEW 15 cm	4011-01
Large	L	3.3 mm	1.7 mm	1 unit per pack	1.7 mm	3.3 mm	35	8 cm	4009-07
Large	L	1.7 mm	1.7 mm	1 unit per pack	1.7 mm	3.3 mm	35	NEW 15 cm	4011-02
Extra Large	XL	4.93 mm	2.5 mm	1 unit per pack	2.5 mm	4.9 mm	25	8 cm	4009-08
Extra Large	XL	2.5 mm	2.5 mm	1 unit per pack	2.5 mm	4.9 mm	25	NEW 15 cm	4011-03



AnastoClip Reusable Atraumatic Forceps 18 cm



AnastoClip Universal Clip Remover 6 units per pack 10 cm  
AnastoClip Universal Clip Remover 6 units per pack 15 cm



# Chevalier Valvulotome- LeMaitre USA

Fast and Effective Valve Removal



**Chevalier Valvulotome** - A simple, yet effective device for the removal of valves during procedures such as in situ bypass. The easy-to-use valvulotome has specially positioned teeth that provide efficient tissue removal, reducing the risk of trauma or damage to the vein wall.

# Eze-Sit<sup>®</sup> Valvulotome-

Flexible, Safe, and Effective



## NEW 1.5mm HYDRO LeMaitre<sup>®</sup> Valvulotome



**Medishare Inc,**

S-439, G.K.1, New Delhi ,



+91-11-29231218

Email-medishare@gmail.com

# VascuTape<sup>®</sup> Radiopaque Tape Glow 'N Tell Radiopaque Precision In Vascular Imaging

Vascutape is a flexible, medical-grade tape for use during any vascular procedure involving fluoroscopy or radiography. Providing a simple way to cross-reference precisely between the inside and outside of a patient's body and to accurately size or locate tributaries or lesions beneath the skin, VascuTape Radiopaque Tape can be used for stent placements, identifying fistula locations, PTA procedures during in situ bypass grafting, vena cava filter placement, atherectomy, and balloon catheter insertions.

VascuTape Radiopaque Tape can be applied to the skin before any radiological procedure. The radiopaque markings are visible both to the eye and to an x-ray machine or fluoroscope, allowing for cross-referencing precisely between the inside and outside of a patient's body. Different types of VascuTape are available for different procedures. Glow 'N Tell<sup>®</sup> (30 cm), LeMaitre's original tape product, has centimetre markings and the LeMaitre<sup>®</sup> Stent Guide (270 mm) has millimetre markings.



Reduces the need for exposure to contrast media.

Flexible tape follows patient contours and does not tear under normal use.

Smaller skin incisions; less wound necrosis. Applies directly to patient's skin with more accurate lesion locating, more precise stent and catheter sizing, and reduced contrast injections. Will move with the patient.

## Application : -

- \* Percutaneous Transluminal
- \* Angioplasty (PTA):
- \* Stent placement
- \* In Situ Vein Bypass
- \* Plaque Excision




# Aortic Occlusion Catheters

## LeMaitre USA

Temporary Atraumatic Vessel Occlusion and Safe, Rapid Control of Bleeding


- Pruitt Peripheral Occlusion Catheters - 3F, 4F, 5F
- Pruitt Irrigation/Occlusion Catheter - 9F
- Pruitt Aortic Occlusion Catheter - 12F
- LeMaitre Aortic Occlusion Catheters - 8F



# Distal Perfusion Catheter

## LeMaitre USA


For Organ Perfusion During Open Aortic Repair Procedures





# SILICON VESSEL TIES / VESSEL LOOPS

## BATRIK MEDICAL



**BLUE**

**RED**

**WHITE**

**YELLOW**



**Size**

**Mini ( 0.2cm \* 40cm)**

**Maxi (0.3 cm \* 40cm)**

Vessel loops are used in occluding, retracting and identifying arteries, veins, tendons and nerves in surgical procedures. Our vessels loops are extruded from medical grade silicone, with the exception of the Richard – Allan vessel loops, which are made of high performance, medical grade elastomeric. All vessels loops are latex free and radiopaque. Size Mini ( 0.2cm \* 40 cm & Maxi (0.3 cm \* 40cm)

**Cardiology**

# Cardiology

MEDISHARE USA, INC.

## NIKO RANGE OF ECG ELECTRODES


- NIKO - MS 005- Eye ball round electrodes(Adult) Liquid gel
- NIKO - MS 035- Eye ball square electrodes Liquid gel
- NIKO - MS 010- Round electrodes (Pead) Solid gel



slide no.3



## Inflation Device




## AngioSafe Guide wires

### Defibrillation Electrode

OBS Defibrillation Pads are produced by high quality of hydrogel with excellent adhesivity through our unique automatical production line. It's hydrogel is overlapping design which can avoid the exposure of conductive film and provide homogeneous and uniform distribution of current on the surface to avoid burn happen.

OBS Defibrillation Pads are equipped with different connectors compatible with different brand of defibrillators.

- 1 Defibrillation Pad Model: **OBS-DE/P**  
The plug Model: 303A5030  
Connector compatible with: Philips Heartstart FR2, FR2 (+), FR3, Heartstream Forerunner, Heartstart MRx, Heartstart XL, Heartstart XL+
- 2 Defibrillation Pad Model: **OBS-DE/M**  
The plug Model: 303A5035  
Connector compatible with: Medtronic Physio-Control, Osato Baxan, Cardioline, Mindray
- 3 Defibrillation Pad Model: **OBS-DE/W**  
The plug Model: 303A3025  
Connector compatible with: Nihon Kohden  
With pad adaptor JC-755V, TEC-5500 Series, TEC-7600 Series, TEC-7700 Series  
With pad adaptor JC-855V, TEC-8300 Series
- 4 Defibrillation Pad Model: **OBS-DE/Z1**  
The plug Model: 303A2545  
Connector compatible with: Zoll Medical Corp.
- 5 Defibrillation Pad Model: **OBS-DE/Z2**  
The plug Model: 303A6030  
Connector compatible with: Zoll E Series, M Series, R Series, AED PLUS, AED PRO
- 6 The plug Model: 303A1525  
Connector compatible with: CU Medical Systems, Inc., Cmos Drake, Propetti
- 7 The plug Model: 303A2060  
Connector compatible with: Matrax Primedic





# NIKO RANGE OF ECG ELECTRODES

- NIKO Eyeball Shape, liquid gel ECG Electrode – MS005
- NIKO Round Shape, Solid gel ECG Electrode – MS010
- NIKO Square Shape, liquid gel ECG Electrode – MS035
- NIKO Pediatric Electrode – MS004
- NIKO MRI Electrode – MS005 RT

1. I.C.U. / C.C.U.
  2. O.T. /EMERGENCY
  3. HEART STATION
- TMT (Tread Mill)  
Echo  
Holter 24hrs. E.C.G.  
**AND OPERATION THEATER**  
1 Millivolt of Current  
Skin Friendly  
Silver/ Silver Chloride Sensor  
**Meets All ECG Monitoring Requirement**



MS005



MS004



MS035





## NIKO MRI ELECTRODES (CATHLAB ELECTRODE) Niko MRI Electrodes for MRI Centers

For Angiography & Angioplasty Procedure

Carbon Stud

Avoid Shadow Formation

Conducts 1 mv of Current

AGCL Sensor



## NIKO NEO NATAL ELECTRODES These are Hypoallergenic

Meant for children upto three years

Special Softer Adhesive

Reliable & Economical



# ANGIOSAFE

FOR ANGIOGRAPHY AT RELIABLE HAND

The Dispenser, to fit its guide wires



## Angiosafe wires

1. Special PTFE coating
2. Teflon which can glide
3. Two lengths – **150 & 260cm**
4. **Two Shapes – J tip & Straight**
5. **Various diameters - (0.025,0.032,0.035,0.038)**
6. CE marked with special dispenser
7. Costs 60% less than any other wires
8. Angiosafe is safe for all Angio procedures

## HYDROPHILLIC NITINOL GUIDEWIRE



### *Special product with double ended flexible tip*


*One end angled standard tip, other  
end straight soft tip. Stiff shaft  
stiffness. Available in two lengths:  
150 and 180 cm.*

**Available 150 ,  
180 , 220 , 260 J &  
St Tip**



# Defibrillation Electrode

303A5030		Connector compatible with: Phillips Heartstart: FR2, FR2 (+), FR3, Heartstream Forerunner, Heartstart MRx, Heartstart XL, Heartstart XL+
303A5035		Connector compatible with: Medtronic Physiocontrol, Osatu Bexen, Cardioline, Mindray, JOUSING
303A4530		Connector compatible with: Nihon Kohden:
303A2545		Connector compatible with: Zoll Medical Corp
303A6030		Connector compatible with: Zoll: E Series, M Series, R Series, AED PLUS, AED PRO
303A1525		Connector compatible with: CU Medical Systems, Inc, Cmos Drake, Progetti
303A2060		Connector compatible with: <b>Metrax Primedic</b>
303A4531		Connector compatible with: Nihon Kohden
303A7035		Connector compatible with: CU Medical Systems, Inc, Philips
303A4525		Connector compatible with: Schiller
303A4035		Connector compatible with: Mediana



# NIKO INFLATION DEVICE – 40 ATM




Disposable balloon inflation devices  
S type 30ml 40atm with stopcock

Specification: S-type, 30ml, 40atm/bar, HP  
Stopcock

## Features:

1. Square shape handle applies to different sizes of hands.
2. Chinese knots on handle prevent slippery on holding.
3. Soft cover on deflation release button reduces deflation vibration.

# ETHOX - Infu-Surg<sup>®</sup> Pressure Infusion Bag



## Specifications of Pressure bags

1. For emergencies we offer special bags 500ml & 1000 ml
2. Disposable & reusable
3. With manometer & without manometer
4. Special material that does not leak





# CATHLAB / X-RAY / RADIATION PROTECTION ACCESSORIES

(DR GOOS SUPREMA- GMBH / AADCO MEDICAL-USA/ XENOLITE USA)



- Aprons are used to protect from Tumor, Cataract, Hyper & Hypothyroidism
- FDA approved aprons
- Zero weight
- Customized aprons made to your shape and size
- Keeps you safe
- Enhance your working hours in lab
- Lead free

