

MASTER BOOK

2023



HEART CARE INDIA CO

No 21/11, 2nd Floor, VC Garden 2nd Street, Mandaveli, Chennai-600 028.

Phone No: 044-24954590/9962514402


Mail Id: heartcare1995@gmail.com

Web: www.heartcareindia.com


Cardiac Vascular Cardiology

CARDIAC SURGERY

CHASE HEART STABILIZER / SVR











CHASE - INTRACORONARY VESSEL SHUNT







Pruitt F3® Balloon






NIKO Flexible Pen



CARDIAC SURGERY












OTW – Embolectomy







Viper II Vacuum Stabilization System



Complete Line to Beating Heart Stabilization Technology

The Viper II Vacuum Stabilization System combines vacuum stabilization with an apical positioner and blower mister visualization.





BLM-700MS Viper II Blower Mister

Viper II Features

- Low Profile Stabilization for Access to Anastomotic Site
- Malleable Footpads provide Optimal Positioning
- Apical Positioning for Access to All Arteries
- Blower Mister Offers Clear View and Access
- Stabilizer and Apical Positioner Fit Other Sternal Retractors



HEART CARE INDIA CO

No 21/11, 2nd Floor, VC Garden 2nd Street,

Mandaveli, Chennai-600 028.


Phone No: 044-24954590/9962514402

Mail Id: heartcare1995@gmail.com
Web: www.heartcareindia.com

MEDTRONIC equivalent to OCTOPUS 3 Stabilizer and URCHIN Heart Positioner

STB-700VS: Viper II Vacuum Stabilizer, M base API-700VS: Viper II Apical Lift Positioner, M Base

Designed to fit on the Medtronic or Chase Retractor





Vacuum Foot Pods and Vacuum Apical Positioner are identical on both "M" and "Q" bases. They are more equivalent to the Maquet system.



Disposable Retractor Rails to fit on the Maquet Retractor Sold as part of the STB-800VS: Viper II Vacuum Stabilizer with Q Base

New Codes: MAQUET EQUIVALENT

STB-800VS, Viper II Stabilizer with Q Base, includes Retractor Rails (see next image) API-800VS, Viper II Apical Lift Positioner with Q Base

Equivalent to ACROBAT SUV, MODEL: OM-9000 Stabilizer and EXPOSE 4 Apical Positioner

Fits only on the Maquet Retractor


CHASE INTRACORONARY LONG 30 MM SHUNTS





HEART CARE INDIA CO

No 21/11, 2nd Floor, VC Garden 2nd Street, Mandaveli, Chennai-600 028. Phone No: 044-24954590/9962514402 Mail Id: <u>heartcare1995@gmail.com</u>







- Helps prevent Arrythmias
- · Control arterial bleeding
- Faster & more accurate
- Soft Silicone Tip

BULB Tips

- Minimize Trauma
- Obstruction –Free flow Path
- Improve Blood Flow
- Flexible Shaft





Product	Ref Code
Tapered Tip Shunts 1.00/ 12mm Shaft	SHT-1210S
Tapered Tip Shunts 1.00/ 12mm Shaft	SHT-1213S
Tapered Tip Shunts 1.00/ 12mm Shaft	SHT1215S
Tapered Tip Shunts 1.00/ 12mm Shaft	SHT-1218S
Tapered Tip Shunts 1.00/ 12mm Shaft	SHT-1220S
Tapered Tip Shunts 1.00/ 12mm Shaft	SHT-1225S
Bulb Tip Shunts 1.50/ 12mm Shaft	SHB-1215S
Bulb Tip Shunts 1.75/ 12mm Shaft	SHB-1218S
Bulb Tip Shunts 2.00/ 12mm Shaft	SHB-1220S
Bulb Tip Shunts 2.50/ 12mm Shaft	SHB-1225S
Bulb Tip Shunts 3.00/ 12mm Shaft	SHB-1230S
Bulb Tip Shunts 2.00 / 30mm Shaft	SHB-3020S
Bulb Tip Shunts 2.50 / 30mm Shaft	SHB-3025S
Bulb Tip Shunts 3.00 / 30mm Shaft	SHB-3030S

LeMaitre- Range of Vascular Grafts



LifeSpan® ePTFE Vascular Grafts

- Superb suture retention2
- High burst strength2
- for Haemodialysis



AlboGraft® Polyester Vascular Graft

- Knitted and woven double velour grafts
- Approved collagen impregnation process
- Comfortable sealing
- Excellent suture ability
- Softness and conformability
- Removable external spiral for peripheral grafts



Omniflow® II Biosynthetic Vascular Prosthesis

• Made Up of collagen and a polyester mesh For Hemodialysis where vein is not available

XenoSure® Biologic Patch

Exceptionally Strong, Uniform

LeMaitre'

- · Vascular & Cardiac Repair
- · Suture Line Reinforcement
- · Soft Tissue Deficiency Repair
- NEW Dura Repair

and Easy to Suture

STRONG, UNIFORM AND EASY TO SUTURE

- Biocompatible
- · Excellent tensile and suture retention strength
- · Does not require special sutures
- · Uniform collagen thickness* for easier suturing

450K+ XenoSure Patches Implanted over 5 Years**



^{**} Based on LeMatire Patch units sold in the last 5 years



HEART CARE INDIA CO

No 21/11, 2nd Floor, VC Garden 2nd Street, Mandaveli, Chennai-600 028.


Phone No: 044-24954590/9962514402
Mail Id: heartcare1995@gmail.com
Web: www.heartcareindia.com

^{*} Data on file at LeMatire Vascular, Inc.


The SVR System Left Ventricular Repair System



The SVR, procedure is a comprehensive surgical therapy for patients suffering from ischemic heart failure. Studies have shown that surgically restoring a dilated left ventricle to a more normal elliptical anatomy using the SVR system results in significantly improved ejection fraction and less risk of post-operative mitral regurgitation



Spherical Ventricle Post SVR Undesirable surgical outcome



An SVR technique should result in a normal elliptical ventricle every time.

Key Benefits:

- Standardizes the SVR procedure
- Achieves optimal ventricular function and volume
- Defines the new apex
- Reorients the mitral apparatus and papillary muscles properly

Standardizing Surgical Ventricular Restoration (SVR) for consistent and elliptical results every time

The SVR System includes:

The Mannequin™

Defines the new apex

Acts as a template to create an elliptical ventricle with optimal volume every time Offered in various sizes based on the patient's BSA

Patch Seizer

- •Reduces potential for incorrectly trimming the patch
- •Decreases the procedure time

Ordering Information

Product	Code	Description	Jnits
	SVC-0080S	SVR System with 80cc Mannequin, without patch	1
1	SVC-0085S	SVR System with 85cc Mannequin, without patch	1
	SVC-0090S	SVR System with 90cc Mannequin, without patch	1
	SVC-0100S	SVR System with 100cc Mannequin, without patch	1
	SVC-0110S	SVR System with 110cc Mannequin, without patch	1
111000	SVC-0120S	SVR System with 120cc Mannequin, without patch	1
· 1	SVC-0130S	SVR System with 130cc Mannequin, without patch	1
6 12 1	SVC-0140S	SVR System with 140cc Mannequin, without patch	1
11 6	SVC-0150S	SVR System with 150cc Mannequin, without patch	1
10 10	* Patents US668177	3 and additional patents pending	






CHASE -TRIPOD™

CHASE

Beating Heart Stabilization System



Easy to Use and Clinically Flexible:

The TRIPOD™ Beating Heart Stabilization System combines ease of use with a rock solid, metal articulating arm and ultra-low profile foot pads. The strength and flexibility of the metal articulating arm delivers superior access and stabilization even to the most challenging Anastomotic locations.


TRIPOD Features:

Ease of use
Superior arm strength and stability
Low profile design
Optional Atraumatic vessel occlusion



Pruitt F3[®] Carotid Shunts

Cerebral Perfusion with Dual Balloon Occlusion



Pruitt F3-S Polyurethane Outlying Carotid Shunt

Description	Diameter	Length	Model #
Shunt with T-Port	8F	31 cm	2015-10
Shunt with T-Port	9F	31 cm	2014-10

Soft, flexible polyurethane with atraumatic tips to minimize intimal trauma.

Increased Flow: Average flow is 39% - 104%* higher than leading competitors.

Improved Markings: Centimetres markings span the length of the shunt with a bar denoting centre point for easier identification of placement in carotid artery. Tether: Aids in removal of shunt.

*Flow comparison testing conducted at body temperature at constant pressure..



Flexcel[™] Carotid Shunt

Increased Flexibility and Flow



Flexcel Carotid Shunts (Inlying)	Length	Model # REF
Single Unit Packing		
8F - packing with 5x1 unit 8F shunt	14.5 cm	2020-15
10F - packing with 5x1 unit 10F shunt	14.5 cm	2020-25
12F - packing with 5x1 unit 12F shunt	14.5 cm	2020-35
14F - packing with 5x1 unit 14F shunt	14.5 cm	2020-45


NIKO AORTIC PUNCHES

Single Rotating Head











Aortic Punches: 1. Punches with Precision 2. Comfortable Grip 3. Smooth and easy to operate 4. Economical 4.4 mm

5. Higher number of Usages

NIKO POLYMER CLIP & MICRO CLIPS WITH OPEN CLIP APPLICATOR

Chamara	PLTX220 Medium		В	lue	
A THE LAND OF THE PARTY OF THE	PLTX230 Medium-Large 0		Gr	een	
dillillilli	PLTX240	PLTX240 Large		Pu	rple
Section.	PLTX250 Extra-Large		Gold		
Tides Tides	TCTL1201-24		CardioSet-ST		Red
Titte Titte	TCTL2200-	24	CardioSet-	МТ	Blue



SOFT TISSUE RETRACTOR

Minimally Invasive Surgery

CARDIAC SURGERY / CTVS SURGERY / LAPROSCOPIC-GENERAL SURGERY
GYNAECOLOGY / C-SECTION AND BREAST SURGERY



- 1. Assisted thoracic small incision surgery
- 2. Pediatric Congenital heart repair surgery
- 3. Minimally Invasive Cardio Surgery, valve replacement
- 4. The small incision robot arm collaborative surgery



HEART CARE INDIA CO

No 21/11, 2nd Floor, VC Garden 2nd Street,

Mandaveli, Chennai-600 028.

Phone No: 044-24954590/9962514402

Mail Id: heartcare1995@gmail.com
Web: www.heartcareindia.com

BULLDOG CLAMP





Vascu®Stop



Multi Clamping Pressures


BULLDOG CLAMPS

Prdoduct Code	Products Description/ürün	barcode no	Pressure Gr.	In Box
BMXC-006	Maxi Curved	8680041680106	125/175/200 gr	25
BMX5-007	Maxi Streight	8680041680113	125/175/200 gr	25
BSFC- 008	Standard Curved	8680041682766	125/175/200 gr	25
BSF5-009	Standard Streight	8680041682773	125/175/200 gr	25
BMNC-010	Mini Curved	8680041682452	125/150 gr	25
BMN5-011	Mini Streight	8680041682759	125/150 gr	25
B5TC-029	Std.Soft Coated Curved	8680041680632	175/250 gr	25
B5T5-030	Std.Soft Coated Straight	8680041680625	175/250 gr	25
	7			





SILICON VESSEL TIES / VESSEL LOOPS BATRIK MEDICAL



Size Mini (0.2cm * 40cm) Maxi (0.3 cm * 40cm)

Medishare Inc.

S-439, G.K.1, New Delhi, +91-11-29231218
Email-medishare@gmail.com

Vessel loops are used in occluding, retracting and identifying arteries, veins, tendons and nerves in surgical procedures. Our vessels loops are extruded from medical grade silicone, with the exception of the Richard – Allan vessel loops, which are made of high performance, medical grade elastomeric. All vessels loops are latex free and radiopaque. Size Mini (0.2cm * 40 cm & Maxi (0.3 cm * 40cm)

The New Innovative ((70mm-200mm)

NIKO Flexible Cautery with LED/USB

Indications:

- 1. Onco Surgery / Abdominal Whipple
- 2. Liver Cancer Colorectal
- 3. Upper Spine
- 4. Hip Replacement / Hip Resection
- 5. Pelvic Resection / Hysterectomy
- 6. MICS, CABG, Neuro Endoscopy, LIMA
- 7. Kidney Transplant Surgery

Features:

- I. Easy to Operate
- 2. USB Connection with Long Cable
- 3. Ergonomically Design
- 4. Light in Weight
- 5. Adjust to any length
- 6. Long Blade (70mm-200mm)
- 7. Power Source 220ac/ Power Banks/ Laptop or any USB outlet
- 8. On the Spot Illumination

Flexible Pen



Flexible Pen with light



Defibrillation Electrode


303A5030	Connector compatible with: Philips Heartstart: FR2, FR2 (+), FR3, Heartstream Forerunner, Heartstart MRx, Heartstart XL, Heartstart XL+
303A5035	Connector compatible with: Medtronic Physiocontrol, Osatu Bexen, Cardioline, Mindray, Jousing
303A4530	Connector compatible with: Nihon Kohden:
303A2545	Connector compatible with: Zoll Medical Corp
303A6030	Connector compatible with: Zoll: E Series, M Series, R Series, AED PLUS, AED PRO
303A1525	Connector compatible with: CU Medical Systems,inc , Cmos Drake, Progetti
303A2060	Connector compatible with: Metrax Primedic
303A4531	Connector compatible with: Nihon Kohden
303A7035	Connector compatible with: CU Medical Systems ,inc,Philips
303A4525	Connector compatible with: Schiller
303A4035	Connector compatible with: Mediana

Medishare Inc,

S-439, G.K.1, New Delhi, medishare@gmail.com









Niko

Electrosurgical Cautery Pad / Diathermy Pad /With & without **Cables Valleylab / Bertchold / Erbe Connectors**



STANDARD & REM



BUTTERFLY SHAPE




INFANT & NEONATAL



CONNECTORS







WITHOUT CABLE PAD


- Plates are available in two types SPLIT & STANDARD, These are hydro gel combined, very cool & comfortable
- These are maximum temperature 4°C have no burning or marks on the skin
- 3. These are adaptable resistance.
- These are available for Infant, Pediatric & Adult, with & without cable
- These are most economical & use to avoid shock to the patient during surgery.

NIKO Range of Electrosurgical Pencils

NIKO PEN



NIKO PEN WITH NON STICK BLADE



NIKO PRO PENCIL WITH CERAMIC BLADE




NIKO REUSABLE PENCIL



NIKO FLEXIBLE PEN



NIKO FLEXIBLE WITH LIGHT



NIKO SMOKE PENCIL



NIKO PENCIL WITH TIP CLEANER




Vascular Surgery



Pruitt F3[®] Carotid Shunts

Cerebral Perfusion with Dual Balloon Occlusion



Pruitt F3-S Polyurethane Outlying Carotid Shunt

Description	Diameter	Length	Model #
Shunt with T-Port	8F	31 cm	2015-10
Shunt with T-Port	9F	31 cm	2014-10

Soft, flexible polyurethane with atraumatic tips to minimize intimal trauma.

Increased Flow: Average flow is 39% - 104%* higher than leading competitors.

Improved Markings: Centimetres markings span the length of the shunt with a bar denoting centre point for easier identification of placement in carotid artery. Tether: Aids in removal of shunt.

*Flow comparison testing conducted at body temperature at constant pressure..




Flexcel[™] Carotid Shunt

Increased Flexibility and Flow



Length	Model # REF
14.5 cm	2020-15
14.5 cm	2020-25
14.5 cm	2020-35
14.5 cm	2020-45
	14.5 cm 14.5 cm 14.5 cm

LeMaitre® Clot Management Catheters



LeMaitre Embolectomy Catheters -

- Concentric balloons with progressive Inflation
- Radiopaque catheter shaft
- French size colour coding Depth markings each 10 cm
- Robust Navigation

Over-the -wire Latex **Embolectomy Catheter**

Single Lumen Latex Embolectomy Catheter

Novasil Single Lumen Catheter

Latis Graft Cleaning Catheter

Unique Latis mesh allows for efficient cleaning and effective removal of clots from synthetic grafts. The traction balloon engages tough adherent thrombi and feedback for maintains tactile effective treatment



Shaft design ensures stability in center of vessel to minimize potential damage to the vessel intima. Non-fragmenting balloon enhances clinical outcomes.



XenoSure® Biologic Patch

Exceptionally Strong, Uniform and Easy to Suture

LeMaitre'

- · Vascular & Cardiac Repair
- · Suture Line Reinforcement
- · Soft Tissue Deficiency Repair
- NEW Dura Repair

STRONG, UNIFORM AND EASY TO SUTURE

- Biocompatible
- · Excellent tensile and suture retention strength
- · Does not require special sutures
- · Uniform collagen thickness* for easier suturing

450K+ XenoSure Patches Implanted over 5 Years**



^{**} Based on LeMatire Patch units sold in the last 5 years



HEART CARE INDIA CO

No 21/11, 2nd Floor, VC Garden 2nd Street, Mandaveli, Chennai-600 028. Phone No: 044-24954590/9962514402 Mail Id: heartcare1995@gmail.com

Web: www.heartcareindia.com

^{*} Data on file at LeMatire Vascular, Inc.

LeMaitre- Range of Vascular Grafts



LifeSpan® ePTFE Vascular Grafts

- Superb suture retention2
- High burst strength2

for Haemodialysis



AlboGraft® Polyester Vascular Graft

- Knitted and woven double velour grafts
- Approved collagen impregnation process
- Comfortable sealing
- Excellent suture ability
- Softness and conformability
- Removable external spiral for peripheral grafts



Omniflow® II Biosynthetic Vascular Prosthesis

• Made Up of collagen and a polyester mesh For Hemodialysis where vein is not available

VenaBLOCK -Varicose Vein Closure System

VenaBLOCK® offers a unique treatment for venous insufficiency. Treatment is simply delivering cyanoacrylate inside of the dysfunctioning vein segment. (VSM, VSP, or perforator veins) continuously while applying pressure over target vein segments



Catheter Length: 100 cm (PTFE)

Catheter Material: Marked on every 2 cm

OD: 2mm

Catheter Diameter 15 cm

Catheter Shaft: Braided, laser Guided

Marker: Double Lumen, ID 0.5mm for agent,

ID 0.9mm for laser

TemREN - Atherectomy System

Ready to use set content provides a practical application option. Rotation speed can be adjusted easily and provide safe space to surgeon.

Consistent Solution For Peripheral Arterial Disease


Ready to use set content provides a practical application option. No additional installation equipment required before the process.


Rotation speed can be adjusted easily and provide safe space to surgeon.





AnastoClip GC For Interrupted Vascular Anastomosis - AV Fistulas





Key Benefits

- Interrupted, compliant anastomoticline
- Securely grips into vessel walls
- Fewer revisions to maintain functional patency1
- Reduced procedure time2

Applications:

- AV fistulas
- Peripheral bypass vein grafting
- Arteriotomy and venotomy
- Vessel repair

		2,3 mm						
Medium	M		1 unit per pack	1.1 mm	2.3 mm	35	8 cm	4009-06
Medium	M	1.1 mm 3.3 mm	1 unit per pack	1.1 mm	2.3 mm	35	NEW 15 cm	4011-01
Large	L		1 unit per pack	1.7 mm	3.3 mm	35	8 cm	4009-07
Large	L	1.7 mm 4.93 mm	1 unit per pack	1.7 mm	3.3 mm	35	NEW 15 cm	4011-02
Extra Large	XL		1 unit per pack	2.5 mm	4.9 mm	25	8 cm	4009-08
Extra Large	XL	2.5 mm	1 unit per pack	2.5 mm	4.9 mm	25	NEW 15 cm	4011-03



Chevalier Valvulotome- LeMaitre USA

Fast and Effective Valve Removal

Eze-Sit® Valvulotome-

Flexible, Safe, and Effective



Chevalier Valvulotome - A simple, yet effective device for the removal of valves during procedures such as in situ bypass. The easy-to-use valvulotome has specially positioned teeth that provide efficient tissue removal, reducing the risk of trauma or damage to the vein wall.

Medishare Inc,

S-439, G.K.1, New Delhi, +91-11-29231218
Email-medishare@gmail.com






VascuTape® Radiopaque Tape Glow 'N Tell Radiopaque

Precision In Vascular Imaging

Vascutape is a flexible, medical-grade tape for use during any vascular procedure involving fluoroscopy or radiography. Providing a simple way to cross-reference precisely between the inside and outside of a patient's body and to accurately size or locate tributaries or lesions beneath the skin, VascuTape Radiopaque Tape can be used for stent placements, identifying fistula locations, PTA procedures during in situ bypass grafting, vena cava filter placement, atherectomy, and balloon catheter insertions.

VascuTape Radiopaque Tape can be applied to the skin before any radiological procedure. The radiopaque markings are visible both to the eye and to an x-ray machine or fluoroscope, allowing for cross-referencing precisely between the inside and outside of a patient's body. Different types of VascuTape are available for different procedures. Glow 'N Tell® (30 cm), LeMaitre's original tape product, has centimetre markings and the LeMaitre® Stent Guide (270 mm) has

millimetre markings.



Reduces the need for exposure to contrast media.

Flexible tape follows patient contours and does not tear under normal use.

Smaller skin incisions; less wound necrosis. Applies directly to patient's skin with more accurate lesion locating, more precise stent and catheter sizing, and reduced contrast injections. Will move with the patient.

Application: -

- *Percutaneous Transluminal
- *Angioplasty (PTA):
- *Stent placement
- *In Situ Vein Bypass
- *Plaque Excision

Aortic Occlusion Catheters LeMaitre USA


Temporary Atraumatic Vessel Occlusion and Safe, Rapid Control of Bleeding

- Pruitt Peripheral Occlusion Catheters 3F, 4F, 5F
- Pruitt Irrigation/Occlusion Catheter 9F
- Pruitt Aortic Occlusion Catheter 12F
- LeMaitre Aortic Occlusion Catheters 8F




Distal Perfusion Catheter LeMaitre USA

For Organ Perfusion During Open Aortic Repair Procedures



SILICON VESSEL TIES / VESSEL LOOPS BATRIK MEDICAL



Size Mini (0.2cm * 40cm) Maxi (0.3 cm * 40cm) Vessel loops are used in occluding, retracting and identifying arteries, veins, tendons and nerves in surgical procedures. Our vessels loops are extruded from medical grade silicone, with the exception of the Richard – Allan vessel loops, which are made of high performance, medical grade elastomeric. All vessels loops are latex free and radiopaque. Size Mini (0.2cm * 40 cm & Maxi (0.3 cm * 40cm)

Cardiology

Cardiology

MEDISHARE USA, INC.

NIKO RANGE OF ECG ELECTRODES

NIKO - MS 005- Eye ball round electrodes(Adult) Liquid gel NIKO - MS 035- Eye ball square electrodes Liquid gel


NIKO - MS 010 - Round electrodes (Pead) Solid gel


























NIKO RANGE OF ECG ELECTRODES

- •NIKO Eyeball Shape, liquid gel ECG Electrode MS005
- •NIKO Round Shape, Solid gel ECG Electrode MS010
- •NIKO Square Shape, liquid gel ECG Electrode MS035
- NIKO Pediatric Electrode MS004
- NIKO MRI Electrode MS005 RT







MS005 MS004 MS035

- 1. I.C.U. / C.C.U.
- 2. O.T. /EMERGENCY
- 3. HEART STATION

TMT (Tread Mill)

Echo

Holter 24hrs. E.C.G.

AND OPERATION THEATER

1 Millivolt of Current Skin Friendly

Silver/ Silver Chloride Sensor

Meets All ECG Monitoring Requirement







NIKO MRI ELECTRODES (CATHLAB ELECTRODE) Niko MRI Electrodes for MRI Centers

For Angiography & Angioplasty Procedure

Carbon Stud

Avoid Shadow Formation

Conducts 1 mv of Current

AGCL Sensor

NIKO NEO NATAL ELECTRODES

These are Hypoallergenic

Meant for children upto three years


Special Softer Adhesive

Reliable & Economical

ANGIOSAFE

FOR ANGIOGRAPHY AT RELIABLE HAND

The Dispenser, to fit its guide wires



Angiosafe wires

- 1. Special PTFE coating
- 2. Teflon which can glide
- 3. Two lengths **150 & 260cm**
- 4. Two Shapes J tip & Straight
- 5. Various diameters (0.025,0.032,0.035,0.038)
- 6. CE marked with special dispenser
- 7. Costs 60% less than any other wires
- 8. Angiosafe is safe for all Angio procedures

HYDROPHILLC NITINOL GUIDEWIRE



Special product with double ended flexible tip

One end angled standard tip, other end straight soft tip. Stiff shaft stiffness. Available in two lengths: 150 and 180 cm.

> Available 150, 180, 220, 260 J & St Tip

Defibrillation Electrode


303A5030	Connector compatible with: Philips Heartstart: FR2, FR2 (+), FR3, Heartstream Forerunner, Heartstart MRx, Heartstart XL, Heartstart XL+
303A5035	Connector compatible with: Medtronic Physiocontrol, Osatu Bexen, Cardioline, Mindray, Jousing
303A4530	Connector compatible with: Nihon Kohden:
303A2545	Connector compatible with: Zoll Medical Corp
303A6030	Connector compatible with: Zoll: E Series, M Series, R Series, AED PLUS, AED PRO
303A1525	Connector compatible with: CU Medical Systems,inc , Cmos Drake, Progetti
303A2060	Connector compatible with: Metrax Primedic
303A4531	Connector compatible with: Nihon Kohden
303A7035	Connector compatible with: CU Medical Systems ,inc,Philips
303A4525	Connector compatible with: Schiller
303A4035	Connector compatible with: Mediana

Medishare Inc,

S-439, G.K.1, New Delhi, medishare@gmail.com









NIKO INFLATION DEVICE – 40 ATM



Disposable balloon inflation devices S type 30ml 40atm with stopcock

Specification: S-type, 30ml, 40atm/bar, HP Stopcock

Features:

- 1. Square shape handle applies to different sizes of hands.
- 2.Chinese knots on handle prevent slippery on holding.
- 3. Soft cover on deflation release button reduces deflation vibration.

ETHOX - Infu-Surg® Pressure Infusion Bag



Specifications of Pressure bags

- 1. For emergencies we offer special bags 500ml & 1000 ml
- 2. Disposable & reusable
- 3. With manometer & without manometer
- 4. Special material that does not leak







CATHLAB / X –RAY / RADIATION PROTECTION ACCESSORIES

(DR GOOS SUPREMA- GMBH / AADCO MEDICAL-USA/ XENOLITE USA)





- Aprons are used to protect from Tumor, Cataract, Hyper & Hypothyroidism
- FDA approved aprons
- Zero weight
- Customized aprons made to your shape and size
- Keeps you safe
- Enhance your working hours in lab
- Lead free





